



**UFI™ ULTRASONIC INSTRUMENTS**

*have been designed for multiple clinical applications and are extremely useful in the access cavity. Whether cleaning, finishing or refining the access cavity, flaring the orifice or smoothing and finishing various restoratives, you're sure to find our instruments indispensable. Our double diamond coating provides superb cutting efficiency and is resistant to Sodium Hypochlorite. (NaCOI)*



*The UFI™ Ultrasonic Advantage is Clear:*

**ACCESS  
VISION  
CUTTING EFFICIENCY  
MULTIPLE TIP CONFIGURATIONS  
WET OR DRY  
SUPERB CRAFTSMANSHIP  
DOUBLE DIAMONDS**




**SAN DIEGO  
SWISS  
MACHINING  
INCORPORATED**

*Universal Finishing Instruments with ultrasonic energy is sure to bring your practice some...*

*Good Vibrations!*



## UFI DOUBLE DIAMOND COATING



*The UFI™ instruments fit most major brands of piezo electric ultrasonic generators having an M3 x 0.6mm thread.*

**\$65.00 PER TIP**

**\$345.00 6 - TIP SET**