

The OUTPOST™

POST REMOVAL SYSTEM


The OUTPOST™ Removal System was designed by San Diego Swiss Machining to get out of difficult situations without damaging the root and tooth.

The OUTPOST™ Removal System features 5 trephine core drills with corresponding taps to engage and grip the post. Once the tap is engaged, the removal pliers complete the process.

The OUTPOST™ is designed to completely engage and remove even the most intransigent post while preserving dentition.



SD SWISS innovative Endo OUTPOST™ will assist you with a retreat case. We have combined patient and doctor comfort by adding light weight pliers and strong taps that can stand up to the pressure for predictable results, we believe our OUTPOST™ Removal System is the best on the market today.



OUTPOST™ Removal
Trephine and Tap # 1



OUTPOST™ Removal
Trephine and Tap #2



OUTPOST™ Removal
Trephine and Tap #3



OUTPOST™ Removal
Trephine and Tap #4



OUTPOST™ Removal
Trephine and Tap #5

SAN DIEGO
 **SWISS**

To Order Toll Free:
(855) 737-9477 • (619) 470-2375
To Order Online: www.sdswiss.com
info@sdswiss.com • We Ship Worldwide

Patent Pending